

# Intracellular Cytokine FACS-Staining Protocol

Prepared by D.K. Mitra, S.C. De Rosa, N. Watanabe

Updated June 18, 1997

## Stimulation

1. Separate PBMC using Ficol-Paque and wash 2.5 times with 10 ml culture medium
2. For prestaining with CD4 and CD62L, suspend cells (at least  $2 \times 10^6$ , but  $4 \times 10^6$  is better) in 100 $\mu$ l culture medium, add stain, and incubate for 15 min at RT.
  - a) Half the cells will be unstimulated: add 50 $\mu$ l of stained cells to 950 $\mu$ l of culture medium in 24 well plate
  - b) Half the cells will be stimulated: add 1 $\mu$ l of metalloproteinase inhibitor (1:1000=10 $\mu$ M, stock is 10mM) to the 50 $\mu$ l of remaining cells; mix gently, and incubate a few minutes. In the meantime, to 947 $\mu$ l of culture medium in the 24 well plate, add 1 $\mu$ l PMA (diluted 1 $\mu$ l to 1000 from stock), 1 $\mu$ l ionomycin, 1 $\mu$ l monensin, and mix. If you choose to culture instead in 2ml, then double these amounts. Then add the cells, gently mix using P1000 pipette, and incubate at 37°C for 6 hours.
3. For prestaining with  $\gamma\delta$  reagents, the procedure is the same except that after the staining, the cells must be washed once with large volume (3-5ml). Also, metalloproteinase inhibitor is not used. **Special note:** The unstimulated control cultures should not be prestained with the  $\gamma\delta$  reagents; instead the  $\gamma\delta$  reagents should be included with the other surface reagents at the time of staining. (The prestained unstimulated  $\gamma\delta$  cells lose staining over 6 hrs in culture.)
4. For the cells that will be stained for **perforin**, no stimulation or prestaining is required. Set aside an appropriate number of cells and leave them on ice or at room temp during the 6 hr time period so that all cells can be stained simultaneously.
5. Harvest cells into conical tubes and place them on ice at the end of culture period
6. Pellet cells
7. Wash the cells once with ice-cold PBS/BSA/Azide (2ml)



### Surface staining:

- EMA  
10 → 490  $\mu$ l  
1 → 20  $\mu$ l  
total  
1/1000 dil
8. Suspend cells in ice-cold PBS/BSA/Azide (50 $\mu$ l for each test =  $1-2 \times 10^6$  cells)
  9. Transfer cells to 96-well plate containing the surface antibodies (if using EMA to exclude dead cells, include with surface stains, final conc. = 5 $\mu$ g/ml). If using PI (e.g., with perforin), include PI with the surface stains; use 50x, twice the usual concentration.
  10. Stain for 15 min on ice in dark
  11. Add 150 $\mu$ l PBS and spin. Wash once with 200  $\mu$ l PBS and spin.
  12. Resuspend in 100 $\mu$ l PBS (containing no protein). If using EMA, expose to light for 5 to 10 minutes on ice at this point before fixation.

### Fixation:

13. Add 100 $\mu$ l of 4% Formaldehyde and mix well (final conc. = 2%)
14. Let the cells stand at RT in dark for 20 min
15. Pellet the cells
16. Wash the cells twice with ice-cold PBS/BSA/Azide buffer
17. Pellet the cells and resuspend in 150  $\mu$ l of permeabilization buffer. Mix gently with multichannel pipette
18. Incubate cells at room temperature for 10 min
19. Pellet cells (aspirate very carefully, or flick out supernatant)

### Intracellular staining

20. Resuspend in 25  $\mu$ l per well of permeabilization buffer containing intracellular Abs (cytokine, perforin)
21. Incubate cells at RT for 30 min in the dark
22. Wash twice (2.5x) with 150-200  $\mu$ l of permeabilization buffer. (Note that one wash may be sufficient, but more washes may decrease the background).
23. Wash twice with 200  $\mu$ l of PBS/BSA/Azide buffer (no Saponin)
24. Suspend cells in 200  $\mu$ l of PBS/BSA/Azide buffer (no Saponin) and transfer to FACS tubes. For HIV samples, this final resuspension should be in staining media (without serum) containing 0.5% paraformaldehyde.
25. FACS analysis

## Reagents and Solutions

### 1) PMA + Iono

**\*Stock soln.**

**1mg/ml PMA (DMSO)**

**2 mM Iono (DMSO)**

### 2) Monensine and KB 8301 (metalloproteinase inhibitor)

**Stock soln.**

**2 mM Monensine (ethanol)**

**10 mM KB 8301 (DMSO), as sent supplied by PharMingen, dissolve  
0.5mg powder in 120 $\mu$ l DMSO. Store frozen.**

### 3) PBS/BSA/Azide buffer

<b>Mix:</b>	<b>50 ml</b>	<b>10x PBS, pH 7.4</b>
	<b>450 ml</b>	<b>Cell culture grade H<sub>2</sub>O</b>
	<b>(+) 0.5 ml</b>	<b>1 M Azide</b>
	<b>500 ml</b>	

**Layer 2.5 g of BSA on top of liquid mixture**

**Allow BSA to dissolve at RT without stirring**

**Sterile filter the mixture**

**Store at 4°C**

### 4) Permeabilization buffer:

<b>Mix:</b>	<b>5 ml</b>	<b>10% Saponin in PBS</b>
	<b>(+) 95 ml</b>	<b>PBS/BSA/Azide buffer</b>
	<b>100 ml</b>	

### 5) 10% Saponin

<b>Mix:</b>	<b>5 g</b>	<b>Saponin (Sigma)</b>
	<b>50 ml</b>	<b>PBS, pH7.4</b>

**Place at 37°C until the saponin has dissolved completely**

**Sterile filter the mixture (0.22  $\mu$ )**

**Store at 4°C**